

# Plasmid preparation by alkaline lysis

Under alkaline conditions, high molecular weight DNA is selectively denatured while covalently closed circular plasmid DNA remains double-stranded. When the solution is rapidly neutralized, denatured chromosomal DNA forms an insoluble clot, while plasmid DNA remains in solution (Birnboim and Doly, 1979).

This method yields tens of micrograms of highly pure plasmid DNA from mini-preps (a few milliliters of culture), and up to milligrams from a liter-scale culture.

The plasmids are purified over a silica matrix (*Alternative A*) or by alcohol precipitation (*Alternative B*).

## Risk assessment

▷ Wear gloves, safety glasses, lab coat



Reviewed: May 30, 2023

## Procedures

### >> Alkaline lysis

- |  |  |
|--|--|
| <input type="checkbox"/> Buffer P1 (Resuspension buffer) | <input type="checkbox"/> Buffer N3 (Neutralization buffer) |
| <input type="checkbox"/> Buffer P2 (Lysis buffer)        |  |

(1.) Grow bacteria in 1–5 mL LB (Miller) overnight to stationary phase. Collect the cells by centrifugation.

*Hint:* For culture tubes, spin at  $3\,000 \times g$  for 5 min. For 2 mL centrifuge tubes, a brief spin at  $11\,000 \times g$  is sufficient.

(2.) Resuspend the pellet in 150  $\mu\text{L}$  Buffer P1.

*Hint:* Buffers P1, P2, N3, PB, and PE refer to standard silica-column miniprep kit components (e.g., QIAGEN QIAprep). Equivalent buffers from other suppliers can be substituted. Homemade buffer recipes are available elsewhere.

(3.) Add 250  $\mu\text{L}$  Buffer P2. Mix by gently inverting the tube 6–8 times. Incubate at room temperature for no more than 5 min or until the lysate appears clear.

*Critical:* Do not vortex or pipette to avoid shearing chromosomal DNA.

(4.) Add 350  $\mu\text{L}$  Buffer N3. Mix gently by inversion 6–8 times.

*Critical:* Ensure complete neutralization to precipitate protein and chromosomal DNA. If Buffer P2 contains thymolphthalein, mix until the solution is colorless.

(5.) Centrifuge at  $17\,000 \times g$  for 3 min to pellet precipitated material. Repeat if the lysate remains turbid.

(6.) Transfer the clear supernatant to a fresh tube for purification. Avoid any carryover of precipitate.

### A > Purification of plasmid DNA over a silica matrix

- |   |  |
|---|--|
| <input type="checkbox"/> Buffer PB (Binding buffer) | <input type="checkbox"/> 5 mM Tris hydrochloride |
| <input type="checkbox"/> Buffer PE (Wash buffer)    |  |

(1.) Apply up to 700  $\mu\text{L}$  of supernatant to a silica spin column. Centrifuge at  $11\,000 \times g$  for 1 min, or use a vacuum manifold.

*Hint:* Silica membranes typically bind 15–25  $\mu\text{g}$  DNA. They may be reused to purify more of the same DNA.

(2.) *Optional:* Wash the silica matrix with 450  $\mu\text{L}$  Buffer PB.

*Hint:* Recommended if plasmids are isolated from strains containing high endonuclease activity such as HB101 or JM109.

(3.) Wash with 650  $\mu\text{L}$  Buffer PE.

(4.) Dry the silica matrix by centrifugation at  $11\,000 \times g$  for 2 min.

## Plasmid preparation by alkaline lysis

- (5.) Place the spin column into a clean microcentrifuge tube.
- (6.) Apply 25–55  $\mu\text{L}$  5 mM Tris pH 8.0 to the membrane. Incubate for 1 min, then elute the plasmid by centrifugation at  $11\,000 \times g$  for 1 min. ✂

*Hint:* For higher yield or large plasmids, warm elution buffer to  $70\text{ }^\circ\text{C}$ .

- (7.) *Optional:* Repeat elution using the same or half the volume to improve yield. Pool eluates if desired. +

### B > Purification by alcohol precipitation

- |  |                                    |
|--|------------------------------------|
| <input type="checkbox"/> Isopropyl alcohol       | <input type="checkbox"/> 5 mM Tris |
| <input type="checkbox"/> Buffer PE (Wash buffer) |                                    |

- (1.) Transfer up to 700  $\mu\text{L}$  of supernatant to a clean microcentrifuge tube.
- (2.) Add an equal volume of isopropanol to precipitate the DNA. Vortex to mix thoroughly.
- (3.) Centrifuge at  $11\,000 \times g$  for 5 min. Discard the supernatant and briefly spin to remove residual solvent. ✂  
The pellet will appear nearly transparent.
- (4.) Wash with 650  $\mu\text{L}$  Buffer PE or freshly prepared 70% ethanol without disturbing the pellet.
- (5.) Air-dry the pellet for 2–5 min, or dry in a vacuum centrifuge for 1 min.
- (6.) Resuspend in 25–55  $\mu\text{L}$  5 mM Tris pH 8.0 or water.

### Analyses

- Measure DNA concentration and purity by UV absorbance at 260 nm, 280 nm, and 230 nm.

Nucleic acid	A260 = 1.0	A260/A280	A260/A230
dsDNA	50 ng/ $\mu\text{L}$	1.8–1.9	2.0–2.2

*Quality assurance:* Extinction coefficients vary with nucleotide composition. A260/A280 ratios increase with pH and vary by 0.2–0.3 between instruments. RNA contamination elevates the ratio; protein or phenol contamination depresses it. ◇

- Confirm plasmid identity by restriction digestion 📄 SOP 0033 and agarose gel electrophoresis 📄 SOP 0032.

*Hint:* Use one or more enzymes predicted to linearize the plasmid or release a diagnostic fragment.

### Troubleshooting

#### Alkaline lysis

##### In Step 4:

- Chromosomal DNA does not form a white precipitate but remains viscous.
  - Reduce culture volume and repeat the preparation with gentler handling.

##### In Step 6:

- Chromosomal DNA in the eluate
  - Use cultures within 12–16 h post-inoculation. Older cultures release degraded chromosomal DNA. If not processing immediately, store the bacterial pellet at  $-20\text{ }^\circ\text{C}$ .
  - Mix lysate gently after adding Buffer P2. Limit lysis to 5 min to avoid DNA shearing.
  - After adding Buffer N3, mix immediately and gently to allow complete neutralization.
- RNA contamination in eluate
  - Verify RNase A was added to Buffer P1 before use. Supplement fresh RNase if needed.

## Plasmid preparation by alkaline lysis

---

### *Purification of plasmid DNA over a silica matrix*

*In Step 6:*

- Little or no DNA in eluate
  - Ensure elution buffer is low salt and within pH 7.0 and pH 8.5.

### *Purification by alcohol precipitation*

*In Step 3:*

- No visible pellet after isopropanol precipitation
  - Centrifuge for longer (10 min  $\mu$ N t). Mark the tube orientation before spinning to locate it on the expected side.
  - Add 0.1 vol of 3M sodium acetate pH 5.2 before adding isopropanol to improve precipitation of dilute samples.

## List of references

H. Birnboim and J. Doly, *Nucleic Acids Res.* **7**(6), 1513—1523 (1979).

## Change log

2016-03-15 Rahul Patharkar Alkaline lysis with ethanol precipitation  
2023-05-30 Benjamin C. Buchmuller Adaptation as SOP

---

**Open Protocol — Part of the *Lab Protocols* collection (2025) by B. C. Buchmuller and contributors.** This document is made available under the Creative Commons Attribution Share Alike 4.0 International License. To view a copy of this license, visit <https://creativecommons.org/licenses/by-sa/4.0/>.

For research use only. Provided in good faith, without warranty or liability for any use or results. Users are responsible for compliance with local regulations and institutional policies.

Current when printed. Visit <https://benjbuch.github.io/check/> or scan the QR code to check for updates.



336cef9

