

Purification of nucleic acids from agarose gels

Recovery of DNA fragments from agarose gels is necessary when undesired side-products cannot be removed or inactivated ahead of sensitive downstream applications such as cloning or sequencing.

The DNA can be electro-eluted into a trough excised ahead of the band (McDonnell et al., 1977), or more conveniently, the band can be excised directly.


To remove the agarose, fast flush extraction (Grey and Brendel, 1992) (*Alternative A*) or chaotropic salts such as guanidinium thiocyanate (Matitashvili and Zavizion, 1997) (*Alternative B*) are commonly used. The latter, as a silica-based method, being convenient, at the expense of lower recovery and higher cost. Fast flush extraction is thus preferred for large sample sets.

This is a bench card. Full protocol available online.



Procedures


>> Separation and excision of the desired DNA fragment from agarose gels


- (1.) Resolve the desired DNA fragment on a 0.6–1.2% agarose gel in TAE or lithium acetate buffer. Load at least 200 ng per lane. 

Critical: To avoid carry-over of DNA, use fresh running buffer and leave empty lanes between samples. 


Critical: Avoid TBE buffer if possible, since borate inhibits downstream enzymatic reactions. 


- (2.) Visualize the gel using a blue light or long-wave UV source in a darkroom or gel box.



Critical: Always use long-wave UV (> 400 nm) and minimize exposure. Avoid pre-imaging gels under UV to prevent DNA damage and pyrimidine dimer formation. 

- (3.) Excise the desired band with a clean surgical blade. Trim excess agarose to minimize the slice size. 
- (4.) Transfer the gel slice to a clean microcentrifuge tube.

A > Fast flush extraction

- (1.) Place the gel slice into a spin column. Centrifuge at 15 000 × *g* for 1 min. Collect the flow-through. 
- (2.) Rotate the column 180° and spin again at 15 000 × *g* for 1 min.

Quality assurance: Use a blue light flashlight to check for residual fluorescence. If DNA remains, add 50 μL 5 mM Tris pH 8.0 and repeat the spin. Expect a faint greenish glow in the eluate. 

- (3.) Combine all flow-throughs. Discard the column.
- (4.) *Optional:* Further purify by silica spin column or ethanol precipitation  SOP 0021. 


B > Extraction with guanidinium thiocyanate

 R0135 Buffer QG (Extraction buffer), 0.5 mL Isopropyl alcohol

- (1.) Weigh the gel slice.
- (2.) Add 350 μL Buffer QG per 100 mg gel. If the gel slice weighs less, supplement with electrophoresis buffer to 100 mg.
- (3.) Melt the agarose at 40–60 °C. Mix occasionally. Cool to room temperature.
- (4.) If the phenol red indicator turns red, add 10 μL 3 M sodium acetate pH 5.0.

Purification of nucleic acids from agarose gels

Critical: Do not skip this step—binding will fail at neutral or basic pH.

- (5.) Add 150 μ L isopropanol per 100 mg agarose.
- (6.) Bind DNA to silica  SOP0021. Skip Buffer PB.



List of references

- M. Grey and M. Brendel, *Curr. Genet* **22**(1), 83—84 (1992).
- M. McDonnell, M. Simon, and F. Studier, *J. Mol. Biol.* **110**(1), 119—146 (1977).
- E. Matitashvili and B. Zavizion, *Anal. Biochem.* **246**(2), 260—262 (1997).

 Recipe (available online)  Troubleshooting (available online)  Notes (available online)

Lab Protocols — © B. C. Buchmuller. Licensed CC-BY-SA-4.0. Research use; provided as-is, without warranty. Current when printed. Visit <https://benjbuch.github.io/check/?q=b50402c> or scan the QR code to read the full protocol or to check for updates.

